



# Recombinant Bovine Phosphorylase b kinase gamma catalytic chain, testis/liver isoform (PHKG2)

<b>Product Code</b>	CSB-BP639762BO
<b>Abbreviation</b>	PHKG2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2KJ16
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTLDVGPEDL LPDWAAAKEF YQKYDPKDV I GRGVSSV VRR CVHRATGQEF AVKIMEVTAE RLSPEQLEEV REATRRETHI LRQVAGHPHI ITLIDSYESS SFMFLVFDLM RKGELFDYLT EKVALSEKET RSIMRSLLEA VSFLHNNNIV HRDLKPENIL LDDNMQIRLS DFGFSCHLEP GEKLRELCGT PGYLAPEILK CSMDETHPGY GKEVDLWACG VILFTLLAGS PPFWHRQRIL MLRMIMEGQY QFSSPEWDDR SDTVKDLISR LLQVDPVERL TAEQALQHPF FERCEGSQAW NLTPRQRFRV AVWTVLAAGR VALSAHRIRP LTKSALLRDP YALRPVRRLI DNCAFRLYGH WVKKGEQQNR AALFQHRPPG PFPMMGPEEE GDSATIAEDE AMLVLG
<b>Source</b>	Baculovirus
<b>Target Names</b>	PHKG2
<b>Protein Names</b>	Recommended name: Phosphorylase b kinase gamma catalytic chain, testis/liver isoform Short name= PHK-gamma-T EC= 2.7.11.19 Alternative name(s): Phosphorylase kinase subunit gamma-2 Serine/threonine-protein kinase PHKG2 EC=
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.