



Recombinant Bovine Phosphorylase b kinase gamma catalytic chain, testis/liver isoform (PHKG2)

Product Code	CSB-MP639762BO
Abbreviation	PHKG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KJ16
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MTLDVGPEDL LPDWAAAKEF YQKYDPKDV I GRGVSSVRR CVHRATGQEF AVKIMEVTAE RLSPEQLEEV REATRRETHI LRQVAGHPHI ITLIDSYESS SFMFLVFDLM RKGELFDYLT EKVALSEKET RSIMRSLLEA VSFLHNNNIV HRDLKPENIL LDDNMQIRLS DFGFSCHLEP GEKLRELCGT PGYLAPEILK CSMDETHPGY GKEVDLWACG VILFTLLAGS PPFWHRQRIL MLRMIMEGQY QFSSPEWDDR SDTVKDLISR LLQVDPVERL TAEQALQHPF FERCEGSQAW NLTPRQRFRV AVWTVLAAGR VALSAHRIRP LTKSALLRDP YALRPVRRLI DNCAFRLYGH WVKKGEQQNR AALFQHRPPG PFPMMGPEEE GDSATIAEDE AMLVLG
Source	Mammalian cell
Target Names	PHKG2
Protein Names	Recommended name: Phosphorylase b kinase gamma catalytic chain, testis/liver isoform Short name= PHK-gamma-T EC= 2.7.11.19 Alternative name(s): Phosphorylase kinase subunit gamma-2 Serine/threonine-protein kinase PHKG2 EC=
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.