



Recombinant Mouse Ribosomal protein S6 kinase alpha-1 (Rps6ka1)

Product Code	CSB-MP020464MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18653
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPLAQLKEPW PLMELVPLDP ENGQTSGEAA GLQPSKDEAI LKEISITHHV KAGSEKADPS QFELLKVLGQ GSFQKVLVLR KVTRPDSGHL YAMKVLKKAT LKVRDRVRTK MERDILADVN HPFVVKLHYA FQTEGKLYLI LDFLRGGDLF TRLSKEVMFT EEDVKFYLAE LALGLDHLHS LGIYRDLKP ENILLDEEGH IKLTDGFLSK EADHEKKAY SFCGTVEYMA PEVVNRQGHT HSADWWSYGV LMGKDRKETM TLILKAKLGM PQFLSTEAQS LLRALFKRNP ANRLGSGPDG AEEIKRHIFY STIDWNKLYR REIKPPFKPA VAQPDDTFYF DTEFTSRTPR DSPGIPPSAG AHQLFRGFSF VATGLMEDDG KPRTTQAPLH SVVQQLHGKN LVFSDGYVVK ETIGVGSYSV CKRCVHKATN MEYAVKVIDK SKRDPSEEIE ILLRYGQHPN IITLKDVIYDD GKHVYLVTEL MRGGELLDKI LRQKFFSERE ASFVLHTISK TVEYLHSQGV VHRDLKPSNI LYVDESGNPE CLRICDFGFA KQLRAENGLL MTPCYTANFV APEVLKRQGY DEGCDIWSLG ILLYTMLAGY TPFANGPSDT PEEILTRIGS GKFTLSGGNW NTVSETAKDL VSKMLHVDPH QRLTAKQVLQ HPWITQKDKL PQSQLSHQDL QLVKGAMAAT YSALNSSKPT PQLKPIESSI LAQRRVRKLP STTL
Source	Mammalian cell
Target Names	Rps6ka1
Protein Names	Recommended name: Ribosomal protein S6 kinase alpha-1 Short name= S6K- alpha-1 EC= 2.7.11.1 Alternative name(s): 90 kDa ribosomal protein S6 kinase 1 Short name= p90-RSK 1 Short name= p90RSK1 Short name= p90S6K
Expression Region	1-724
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been



characterized.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.