



Recombinant Bovine Four and a half LIM domains protein 2 (FHL2)

Product Code	CSB-MP640664BO
Abbreviation	FHL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KI95
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MTERFDCHHC EDSLFGKRYV LREEQPYCVA CFEALFASTC EECGKLGICD CKDLSYKDRH WHEACFHCSR CRGSLVDKPF AAKEDQLLCT DCYSQEYSSR CQECKKSIMP GTRKMEYKGS SWHETCFICH RCQQPIGTKS FIPKDSENF VPCYERQYAL QCVQCKKPIT TGGVTYREQP WHRECFVCTA CKKPLSGQRF TSRDEFAYCL GCFCDLYAKK CAGCANPISG LGGTKYISFE ERQWHNDCFN CKKCSLSLVG RGFLTERDDI LCPDCGKDI
Source	Mammalian cell
Target Names	FHL2
Protein Names	Recommended name: Four and a half LIM domains protein 2 Short name= FHL-2
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.