



Recombinant Pig 2'-5'-oligoadenylate synthase 1 (OAS1)

Product Code	CSB-EP640913PI-B
Abbreviation	OAS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29599
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	MELRHTPARD LDKFIEDHLL PNTCFRTQVK EAIDIVCRFL KERCFQGTAD PVRVSKVVKG GSSGKGTTLR GRSDADLVVF LTKLTSFEDQ LRRRGEFIQE IRRQLEACQR EQKFKVTFEV QSPRENPRALSFVLSSPQL QQEVEFDVLP AFDALGQWTP GYKPNPEIYV QLIKECKSRG KEGEFSTCFT ELQRDFLRNR PTKLKSLIRL VKHWYQTCKK THGNKLPPQY ALELLTVYAW EQGSRKTDFFS TAQGFQTVLE LVLKHQKLCI FWEAYYDFTN PIVGRCMLQQ LKKPRPVILD PADPTGNVGG GDTHSWQRLA QEARVWLGYP CCKNLDGSLV GAWTMLQKI
Source	E.coli
Target Names	OAS1
Protein Names	Recommended name: 2'-5'-oligoadenylate synthase 1 Short name= (2-5')oligo(A) synthase 1 Short name= 2-5A synthase 1 EC= 2.7.7.- Alternative name(s): p42 OAS
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.