



Recombinant Danio rerio Olfactomedin-like protein 3A (olfml3a)

Product Code	CSB-EP640938DIL-B
Abbreviation	olfml3a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29RB4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	<p> QQQ ALMDYLERRL LAIEDRISLW HEQTSRYASE LREFKQQMVA QLEGLDKSKE ALRSELDTVG TRVDRVEREM DYLETQNGAQ PCVDVDDKLV EQQVTFVKER NKAKYAKLTD CSDMISSIKA MKILKRVGGA KGMWTKDMGS ATGKVYILNG TDDNTVFEFG TVREFTSSQG TSGATTIQLP SAWRGMGHAI YNNHMYYLKQ GEEMKLIKFD LQKNIIVDSA VFPVKNQLPV YSLNPETFID LAVDEEGLWA IYATQENERH ISLAKIDPKT LDIQQMWDTTP CVRENAEAAF VICGTVYVVY NSKLPSRSRI QCVFDVSDMV NNDEAPHVYF PKRYGTHSSL KYSPVEQLLY AWDDGYQILY KLQLKKKLEV </p>
Source	E.coli
Target Names	olfml3a
Protein Names	Recommended name: Olfactomedin-like protein 3A
Expression Region	18-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.