



Recombinant Phoca hispida Interleukin-1 beta (IL1B)

| | |
|--------------------------|---|
| Product Code | CSB-MP641327PAAZ |
| Abbreviation | IL1B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q2HZH0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Pusa hispida (Ringed seal) (Phoca hispida) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | AP LQSLDCKFRD KNQKSLVLYN SYELRALHLN GSSVNQQAVF RMSFLQGDEN SNKIPVALCI KEKNLYLSCV MKDGGKPTLQL EMLDPRVYPK KKMEKRFVFN KTEVKKILEF ESSQFPNWYI STSKAEAMPV FLGNTKGGQD ITDFTMEFSS |
| Source | Mammalian cell |
| Target Names | IL1B |
| Protein Names | Recommended name: Interleukin-1 beta Short name= IL-1 beta |
| Expression Region | 119-270 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |