



# Recombinant Chicken Protein Wnt-3a (WNT3A)

<b>Product Code</b>	CSB-EP641688CH-B
<b>Abbreviation</b>	WNT3A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2LMP1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SY PIWWSLAIGH QYSSLGTQPI LCGSIPGLVP KQLRFCRNYV EIMPSVAEGV KIGIQECQHQ FRGRRWNCTT VNDSLAIFGP VLDKATRESA FVHAIASAGV AFAVTRSCAE GSATICGCDT RHKGSPGEGW KWGGCSEDVE FGSMVSREFA DARENRPDAR SAMNRHNNEA GRTSIIELMH LKCKCHGLSG SCEVKTCWWS QPDFRVIGDY LKDKYDSASE MVVEKHRESR GWWETLRPKY NFFKAPTEKD LVYYENSPNF CEPNPETGSF GTRDRICNVT SHGIDGCDLL CCGRGHNTRT EKRKEKCHCI FHWCCYVRCQ ECIRVYDVHT CK
<b>Source</b>	E.coli
<b>Target Names</b>	WNT3A
<b>Protein Names</b>	Recommended name: Protein Wnt-3a
<b>Expression Region</b>	19-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 96% amino acid identity to mouse Wnt3A protein, and 84% to human WNT3 protein, another WNT gene product. This gene is clustered with WNT14 gene, another family member, in chromosome 1q42 region.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.