



Recombinant Chicken Protein Wnt-3a (WNT3A)

Product Code	CSB-MP641688CH
Abbreviation	WNT3A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2LMP1
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	SY PIWWSLAIGH QYSSLGTQPI LCGSIPGLVP KQLRFCRNYV EIMPSVAEGV KIGIQECQHQ FRGRRWNCTT VNDSLAIFGP VLDKATRESA FVHAIASAGV AFAVTRSCAE GSATICGCDT RHKGSPGEGW KWGGCSEDVE FGSMVSREFA DARENRPDAR SAMNRHNNEA GRTSIIELMH LKCKCHGLSG SCEVKTCWWS QPDFRVIGDY LKDKYDSASE MVVEKHRESR GWWETLRPKY NFFKAPTEKD LVYYENSPNF CEPNPETGSF GTRDRICNVT SHGIDGCDLL CCGRGHNTRT EKRKEKCHCI FHWCCYVRCQ ECIRVYDVHT CK
Source	Mammalian cell
Target Names	WNT3A
Protein Names	Recommended name: Protein Wnt-3a
Expression Region	19-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 96% amino acid identity to mouse Wnt3A protein, and 84% to human WNT3 protein, another WNT gene product. This gene is clustered with WNT14 gene, another family member, in chromosome 1q42 region.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.