



Recombinant Bovine Neuroguidin (NGDN)

Product Code	CSB-MP642529BO
Abbreviation	NGDN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KII6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	AAPEVLESD LPNAVALLKN LQEQVMVTA QVQTLTKKVQ AKAYPTEKGL SLLEVVDQQL LMYLMDLSHL ILDKASGGSL QGHPAVLRLV EIRTVLEKLR PLDQKLKYQI DKLVKTAVTG SLENDPLRF KPHPSNMMSK LSSEDEEEDE AEEGQSGASG KKSQKGTAKK YVPPRLVPVH YDETEAEREK KRLERAKRRA LSSSVIRELK EQYSDAPEEI RDARHPHVTR QSQEDQHRIN YEESMMVRLS VSKREKGRRK RANVMSSQLH SLTHFSDISA LTGGTPHLDE DQNPTKKRKK IPKKGRKKKG FRRRR
Source	Mammalian cell
Target Names	NGDN
Protein Names	Recommended name: Neuroguidin Alternative name(s): EIF4E-binding protein
Expression Region	2-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.