



# Recombinant *Oryza sativa* subsp. japonica Hexokinase-2 (HXK2)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP642550OFG  |
| <b>Abbreviation</b>      | HXK2   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q2KNB9   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. japonica (Rice)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MRKAAAAAVA AAAAVGVALL VRRQLREAKR WGRADAVLRE LEERCAAPPG<br>RLRQVADAMA VEMHAGLASE GGSCLKMIIS YVDALPSGEE KGVFYALDLG<br>GTNFRVLRVQ LGGKEGRVIK QEHEISIPP HLMTGGSNEL FDFIASSLAK<br>FVASEGEDFH LAEGRQRELG FTFSFPVKQT SIASGTLINW TKGFSIDETV<br>GEDVVELTK ALERQGLDMK VTALINDTIG TLAGGRYDDN DVIAAVILGT<br>GTNAAYVERA NAIPKWHDLL PKSGDMVINM EWGNFRSSHL PLTEFDQALD<br>AESLNPGEQV YEKLISGMYL GEIVRRVLLK MAEEASLFGD EVPPKLIKIPF<br>IIRTPYMSMM HCDRSPDLRT VGAKLKDILG VQNTSLKTRR LVVDVCDIVA<br>KRAAHLAAAG IHGILKKLGR DVPNTDKQRT VIAVDGGLYE HYTIFAECVE<br>STLRDMLGED VSSTIVIKLA KDGSGIGAAL LAAHSQYRE AEEL |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | HXK2   |
| <b>Protein Names</b>     | Recommended name: Hexokinase-2 EC= 2.7.1.1 Alternative name(s):<br>Hexokinase-3  |
| <b>Expression Region</b> | 1-494  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.