



Recombinant Human Protein FAM83E (FAM83E)

Product Code	CSB-YP642651HU
Abbreviation	FAM83E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2M2I3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAASQLAALE GVDSGPRVPG ASPGFLYSEG QRLALEALLS KGAEAFQTCV QREELWPFLS ADEVQGLAAA AEDWTVAKQE PSGMAEGATT TDVDAGSLSY WPGQSEQPAP VLRLGWPVDS AWKGITRAQL YTQPPGEGQP PLKELVRLEI QAAHKLVAVV MDVFTDPDLL LDLVDAATTR WVPVYLLDR QQLPAFLELA QQLGVNPWNT ENVDVRVVRG CSFQSRWRRQ VSGTVREKFV LLDGERVISG SYSFTWSDAR LHRGLVTLT GEIVDAFSLE FRTLYAASCP LPPAPPQKPS VIGGLQRGRS PHRVSRRRSV APASPPPPDG PLAHLAACR VSPATPGPAL SDILRSVQRA RTPSGPPARP SRSMWDL SRL SQLSGSSDGD NELKKSWSGSK DTPAKALMRQ RGTGGGPWGE VDSRPPWGGA LPLPPAHLRL YLSPARRRFG GDAFLKLQEP RGVRPDWDAP RAGLGGQP
Source	Yeast
Target Names	FAM83E
Protein Names	Recommended name: Protein FAM83E
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.