



# Recombinant Human Tectonic-1 (TCTN1)

<b>Product Code</b>	CSB-EP642681HU-B
<b>Abbreviation</b>	TCTN1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2MV58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QTDATPAV TTEGLNSTEA ALATFGTFPS TRPPGTPRAP GPSSGPRPTP VTDVAVLCVC DLSPAQCDIN CCCDPDCSSV DFSVFSACSV PVVTGDSQFC SQKAVIYSLN FTANPPQRVF ELVDQINPSI FCIHITNYKP ALSFINPEVP DENNFDTLMK TSDGFTLNAE SYVSFTTKLD IPTAAKYEYG VPLQTSDSFL RFPSSLTSSL CTDNNPAAFL VNQAVKCTRK INLEQCEEIE ALSMAFYSSP EILRVPDSRK KVPITVQSIV IQSLNKTLTR REDTDVLQPT LVNAGHFSLC VNVVLEVKYS LTYTDAGEVT KADLSFVLGT VSSVVVPLQQ KFEIHFLQEN TQPVPLSGNP GYVVGLPLAA GFQPHKGSGI IQTTNRYGQL TILHSTTEQD CLALEGVRTP VLFGYTMQSG CKLRLTGALP CQLVAQKVKS LLWGQGFPDY VAPFGNSQAQ DMLDWVPIHF ITQSFNRKDS CQLPGALVIE VKWTKYGSLL NPQAKIVNVT ANLISSSFPE ANSGNERTIL ISTAVTFVDV SAPAEAGFRA PPAINARLPF NFFFPFV
<b>Source</b>	E.coli
<b>Target Names</b>	TCTN1
<b>Protein Names</b>	Recommended name: Tectonic-1
<b>Expression Region</b>	23-587
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the tectonic family of secreted and transmembrane proteins. The orthologous gene in mouse is required for formation of most ventral cell types. It functions downstream of smoothed and rab23 to modulate hedgehog signal transduction. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.