



Recombinant *Drosophila pseudoobscura* *pseudoobscura* Drosulfakinins (Dsk)

Product Code	CSB-EP642769DME
Abbreviation	Dsk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q295C6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila pseudoobscura pseudoobscura</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	HQRSTGF
Source	E.coli
Target Names	Dsk
Protein Names	Recommended name: Drosulfakinins Cleaved into the following 3 chains: 1. Drosulfakinin-0 Short name= 2. DSK-0 3. Drosulfakinin-1 Alternative name(s): Drosulfakinin I Short name= DSK-I Drosulfakinin-2 Alternative name(s):
Expression Region	74-80
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.