



Recombinant Bovine Plasma kallikrein (KLKB1), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP643461BO |
| Abbreviation | KLKB1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q2KJ63 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | G CLTQLYHNIF FRGGDVSAMY TPDAQYCQLM CTFHPRCLLF SFLPENSTSD ADKRFGCFLK DSVTGTLP RV SRTGAISGHS LKRCGHQISA CHRSIYKID MRGVNFNASK VRSACECQER CTNNIHCQFF TYATKTFSSA EYRNTCLLKR SPQGTPTRIK VLSDVESGFS LKACGNSKIG CRVDIFQHS A FSDVDVAGII APDAFVCRTI CTYHPSCLFF TFYTNWKT D SQRNVCFLKT SQSGSPSSPT PQENAI SGYS LLTCKQTLPG TEPCHSKIYP QVAFEGEELH VTFVKGVDGC QETCTKMIRC QFFTYSLFPE DCRGEKCKCS LRLSLDGSPT NITYGTQASS GYSLRLCKRG DSRVCTTKRT R |
| Source | Mammalian cell |
| Target Names | KLKB1 |
| Protein Names | Recommended name: Plasma kallikrein EC= 3.4.21.34 Alternative name(s): Fletcher factor Kininogenin Plasma prekallikrein Cleaved into the following 2 chains: 1. Plasma kallikrein heavy chain 2. Plasma kallikrein light chain |
| Expression Region | 20-391 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.