



# Recombinant Bovine [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase, mitochondrial (BCKDK)

<b>Product Code</b>	CSB-BP643467BO
<b>Abbreviation</b>	BCKDK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2KJG8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	STSATDTHHV EMARERSKTV TSFYNQSAID VAAEKPSVRL TPTMMLYSGR SQDGSHLLKS ARYLQQELPV RIAHRIKGFR SLPFIIGCNP TILHVHELYI RAFQKLDFP PIKDQADEAR YCQLVRQLLD DHKDVVTLA EGLRESRKYI EDEKLVRYFL DKTLTSLRGI RMLATHHLAL HEDKPDFVGI ICTRLSPKKI IEKWVDFARR LCEHKYGNAP RVRINGHVAA RFPFIPMPLD YILPELLKNA MRATMESHLD TPYNVPDVVI TIANNIDILV IRISDRGGGI AHKDLDRVMD YHFTTAEAST QDPRISPLFG HLDLHSGGQS GPMHGFGFGL PTSRAYAEYL GGSLRLQSLQ GIGTDVYLRL RHIDGREESF RI
<b>Source</b>	Baculovirus
<b>Target Names</b>	BCKDK
<b>Protein Names</b>	Recommended name: [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase, mitochondrial EC= 2.7.11.4 Alternative name(s): Branched-chain alpha-ketoacid dehydrogenase kinase Short name= BCKD-kinase Short name= BCKDHKIN
<b>Expression Region</b>	31-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.