



Recombinant Bovine [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase, mitochondrial (BCKDK)

Product Code	CSB-EP643467BO-B
Abbreviation	BCKDK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KJG8
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	STSATDTHHV EMARERSKTV TSFYNQSAID VAAEKPSVRL TPTMMLYSGR SQDGSHLLKS ARYLQQELPV RIAHRIKGFR SLPFIIGCNP TILHVHELYI RAFQKLTDFP PIKDQADEAR YCQLVRQLLD DHKDVVTLA EGLRESRKYI EDEKLVRYFL DKTLTSLRGI RMLATHHLAL HEDKPDFVGI ICTRLSPKKI IEKWVDFARR LCEHKYGNAP RVRINGHVAA RFPFIPMPLD YILPELLKNA MRATMESHLD TPYNVPDVVI TIANNIDILV IRISDRGGGI AHKDLDRVMD YHFTTAEAST QDPRISPLFG HLDLHSGGQS GPMHGFGFGL PTSRAYAEYL GGSLRLQSLQ GIGTDVYLRL RHIDGREESF RI
Source	E.coli
Target Names	BCKDK
Protein Names	Recommended name: [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase, mitochondrial EC= 2.7.11.4 Alternative name(s): Branched-chain alpha-ketoacid dehydrogenase kinase Short name= BCKD-kinase Short name= BCKDHKIN
Expression Region	31-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.