



Recombinant Bovine DnaJ homolog subfamily C member 5 (DNAJC5)

Product Code	CSB-MP643681BO
Abbreviation	DNAJC5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29455
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MADQRQRSL S TSGESLYHVL GLDKNATSDD IKKSYRKLAL KYHPDKNPDN PEAADKFKEI NNAHAILTDA TKRNIYDKYG SLGLYVAEQF GEENVNTYFV LSSWWAKALF IFCGLLTCCY CCCCLCCCFN CCCGKCKPKA PEGEETEFYV SPEDLEAQLQ SDEREAADTP IVIQPASATE TTQLTADSHP SYHTDGFN
Source	Mammalian cell
Target Names	DNAJC5
Protein Names	Recommended name: DnaJ homolog subfamily C member 5 Alternative name(s): Cysteine string protein Short name= CSP
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.