



Recombinant Bovine ETS translocation variant 1 (ETV1)

Product Code	CSB-EP644419BO-B
Abbreviation	ETV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KIC2
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MDGFYDQQVP YMVTNNQRGR NCNEKPTNVR KRKFINRDLA HDSEELFQDL SQLQETWLAE AQVPDNDQF VPDYQAESLA FHGLPLKIKK EPHSPCSELG SACSQEQPFK FSYGEKCLYN VSAYDQKPQV GMRPSNPPTP SSTPVSPHHP ASPNSAHTSK PDRAFPALP PSQPIQDSSY PMDHRFRRQL SEPCNSFPPL PTMPREGRPM YQRQMSEPNI PFPPQGFKQE YHDPVYEHNT MVGGAASQSF PPPLMIKQEP RDFAYDSEVP SCHSIYMRQE GFLTHPSRTE GCMFEKGPRQ FYDDTCVVPE KFDGDIKQEP GMYREGPTYQ RRGSLQLWQF LVALLDDPSN SHFIAWTGRG MEFKLEPEE VARRWGIQKN RPAMNYDKLS RSLRYYYEKG IMQKVAGERY VYKFVCDPEA LFSMAFPDNQ RPLLKTDMER HINEEDTVPL SHFDESMAYM PEGGCCNPHP YNEGYYV
Source	E.coli
Target Names	ETV1
Protein Names	Recommended name: ETS translocation variant 1
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.