



# Recombinant Bovine ETS translocation variant 1 (ETV1)

<b>Product Code</b>	CSB-MP644419BO
<b>Abbreviation</b>	ETV1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2KIC2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDGFYDQQVP YMVTNNQRGR NCNEKPTNVR KRKFINRDLA HDSEELFQDL SQLQETWLAE AQVPDNDQF VPDYQAESLA FHGLPLKIKK EPHSPCSELG SACSQEQPFK FSYGEKCLYN VSAYDQKPQV GMRPSNPPTP SSTPVSPHHP ASPNSAHTSK PDRAFPALP PSQPIQDSSY PMDHRFRRQL SEPCNSFPPL PTMPREGRPM YQRQMSEPNI PFPPQGFKQE YHDPVYEHNT MVGGAASQSF PPPLMIKQEP RDFAYDSEVP SCHSIYMRQE GFLTHPSRTE GCMFEKGPRQ FYDDTCVVPE KFDGDIKQEP GMYREGPTYQ RRGSLQLWQF LVALLDDPSN SHFIAWTGRG MEFKLEPEE VARRWGIQKN RPAMNYDKLS RSLRYYYEKG IMQKVAGERY VYKFVCDPEA LFSMAFPDNQ RPLLKTDMER HINEEDTVPL SHFDESMAYM PEGGCCNPHP YNEGYYV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ETV1
<b>Protein Names</b>	Recommended name: ETS translocation variant 1
<b>Expression Region</b>	1-477
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.