



Recombinant Human Exocyst complex component 3-like protein 2 (EXOC3L2)

Product Code	CSB-MP644809HU
Abbreviation	EXOC3L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2M3D2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAALENGELG PLLSPGTLRG LEDECVTVDK AQTRAALLRV LQEDEEHWGS LEDQPSSLAQ DVCELLEEHT ERAPRISQEF GERMAHCCLG GLAEFLQSFQ QRVERFHENP AVREMLPDTY ISKTIALVNC GPPLRALAER LARVGPPSE PAREASASAL DHVTRLCHRV VANLLFQELQ PHFNKLMRRK WLSSPEALDG IVGTLGAQAL ALRRMQDEPY QALVAELHRR ALVEYVRPLL RGRRLCSSAR TRSRVAGRLR EDAAQLQRLF RRLESQASWL DAVVPHLAEV MQLEDTPSIQ VEVGVLVRDY PDIRQKHVAA LLDIRGLRNT AARQEILAVA RDLELSEEGA LSPPRDRAFF ADIPVPRPSF CLSLPLFLGR LPLSRLARPS LACLPRPRPP SLARPRAQR
Source	Mammalian cell
Target Names	EXOC3L2
Protein Names	Recommended name: Exocyst complex component 3-like protein 2 Alternative name(s): HBV X-transactivated gene 7 protein HBV XAg-transactivated protein 7
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.