



Recombinant Bovine Exosome complex component RRP4 (EXOSC2)

Product Code	CSB-MP645613BO
Abbreviation	EXOSC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KID0
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MALEMRLPVA RKPLSESLGR ESKKHLVVPD DTITTDGTFM RGHGTYMGEE KLIASVAGSV ERVKNLICVK ALKTRYNGEV GDIVVGRITE VQQRWQVET NSRLDSVLLL SSMNLPGGEL RRRSAEDELA MRGFLQEGDL ISAEVQAVFS DGAVSLHTRS LKYGKLGQGV LVQVSPSLVK RQKTHFHDLF CGASVILGNN GFIWVYPTPE HKEDDAGGFI ANLEPVSLTD REVISRLRNC IVSLATQRMM LYDTSILYCY EASLPHQIKD ILKPEIMEEI VMETRQRLLLE QEG
Source	Mammalian cell
Target Names	EXOSC2
Protein Names	Recommended name: Exosome complex component RRP4 Alternative name(s): Exosome component 2 Ribosomal RNA-processing protein 4
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.