



Recombinant Solanum tuberosum Patatin-14

Product Code	CSB-EP645766FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2MY37
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	KLEEMVT VLSIDGGGIK GIIPAILEF LEGQLQEVDN NKDARLADYF DVIGGTSTGG LLTAMITTPN ENNRPF ^{AA} AK DIVPFYFEHG PHIFNYSGSI IGPMYDGKYL LQVLQEKLGE TRVHQALTEV AISSFDIKTN KPVI ^{FT} KS ^{NL} AKSP ^{EL} DAKM YDICYSTAAA PIYFPPHYFI THTSNGDIYE FNLVDGGVAT VGDPALLSLS VATRLAQEDP AFSSIKSLDY KQMLLLSLGT GTNSEFDKTY TAQEAAKWGP LRWMLAIQQM TNAASSYMTD YYISTVFQAR HSQNNYLRVQ ENALTGTTTE MDDASEANME LLVQVGETLL KKPVSKDSPE TYEEALKRFA KLLSDRKKLR ASKASY
Source	E.coli
Protein Names	Recommended name: Patatin-14 EC= 3.1.1.-
Expression Region	24-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.