



# Recombinant Bovine Methionine adenosyltransferase 2 subunit beta (MAT2B)

<b>Product Code</b>	CSB-BP645883BO
<b>Abbreviation</b>	MAT2B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q29RI9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VGREKELSI HFVPGDCRLV EEEVNIPNRR VLITGATGLL GRAVYKEFQQ NNWHAVGCGF RRARPKFEQV NLLDSNAVHH IYDFQPHVI VHCAAERRPD VVENHPDAAS QLNVDASGNL AKEAAAIGAF LIYSSDYVF DGTNPPYREE DIPNPLNLYG KTKLEGEKAV LENNLGAAVL RIPVLYGEVE RLEESAVTIM FDKVQFSNKS ANMDHWQQR F PTHVKDVATV CRQLAEKRML DPSIKGTFHW SGNEQMTKYE MACAIADAFN LPSSHLPIT DSPVVGARQP RNAQLDCSRL ETLGIGQRTP FRIGIKESLW PFLIDKRWRQ TVFH
<b>Source</b>	Baculovirus
<b>Target Names</b>	MAT2B
<b>Protein Names</b>	Recommended name: Methionine adenosyltransferase 2 subunit beta Alternative name(s): Methionine adenosyltransferase II beta Short name= MAT II beta
<b>Expression Region</b>	2-334
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein belongs to the methionine adenosyltransferase (MAT) family. MAT catalyzes the biosynthesis of S-adenosylmethionine from methionine and ATP. This protein is the regulatory beta subunit of MAT. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.