



Recombinant Methanospirillum hungatei Enolase 1 (eno1)

Product Code	CSB-EP646030MAAQ-B
Abbreviation	eno1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2FQL9
Product Type	Recombinant Protein
Immunogen Species	Methanospirillum hungatei JF-1 (strain ATCC 27890 / DSM 864 / NBRC 100397 / JF-1)
Purity	>85% (SDS-PAGE)
Sequence	MVQIEKVFGR EIIDSRGNPA VEVDITLVGG YFGRAACPSG ASTGTHEAIE KRDGDFRFFG KGVKKAVEAI NDIRKSLLG MDSTDQPAVD NQLISLDGTD NKGHLGANAI LPVSMAVARA ASQALQVPLY EHLADRTYLL PVPYMNILNG GAHANWQGAD FQEYMIAPVG APDFPEAVRW GCEVYHTLKA ILKKRGLSTG VGDEGGFAPK VSSNRAPLDF ITEAIEMAGY KPGKDISLAL DPASSEFYSD GVYELKSEGK KLSSEEMTGY YEDMVAVYPI ISIEDGLSED DWNGWIAMTK AIGKKVQLVG DDIFVTNTRK ISRGISEKAA NAVLIKLNQI GTVTETISAV TMARNAGWSS MISHRSGETV DSFIADLTVS LGTGQIKTGA PCRGERVEKY NQLMRIHEEL GSRATYAGKK REIRHL
Source	E.coli
Target Names	eno1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.