



Recombinant Bovine DnaJ homolog subfamily A member 2 (DNAJA2)

Product Code	CSB-MP646249BO
Abbreviation	DNAJA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2HJ94
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MANVADTKLY DILGVPPGAS ENELKKAYRK LAKEYHPDKN PNAGDKFKEI SFAYEVLSNP EKRELYDRYG EQGLREGSSG GGGMDDIFSH IFGGGLFSFM GNQSRSRNGR RRGEDMMHPL KVSLEDLYNG KTTKLQLSKN VLC SACSGQG GKSGAVQKCS ACRGRGVRIM IRQLAPGMVQ QM QSVCSDCN GEDEVINEKD RCKKCEGKKV IKEVKILEVH VDKGMKHGQR ITFTGEADQA PGVEPGDIVL LLQEKEHEVF QRDGNDLHMT YKIGLVEALC GFQFTFKHLD GRQIVVKYPP GKVIEPGCVR VVRGEGMPQY RNPFEKGDLY IKFDVQFPEN NWINPDKLSE LEDLLPSRPE VPNIIGDTEE VELQEFDSTR GSGGGQRREA YNDSSDEESS SHHGPGVQC
Source	Mammalian cell
Target Names	DNAJA2
Protein Names	Recommended name: DnaJ homolog subfamily A member 2
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein belongs to the evolutionarily conserved DNAJ/HSP40 family of proteins, which regulate molecular chaperone activity by stimulating ATPase activity. DNAJ proteins may have up to 3 distinct domains: a conserved 70-amino acid J domain, usually at the N terminus; a glycine/phenylalanine (G/F)-rich region; and a cysteine-rich domain containing 4 motifs resembling a zinc finger domain. The product of this gene works as a cochaperone of Hsp70s in protein folding and mitochondrial protein import in vitro.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.