



Recombinant Bovine UBX domain-containing protein 6 (UBXN6)

Product Code	CSB-MP646533BO
Abbreviation	UBXN6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2KIJ6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MKKFFQEIKA DIKFKSAGPG QKLTESVGEK APKEKPSQPP VRQPRQGPTN EAQMAAAAAL ARLEQKQPRA RGPTSQDSIR NQVRKELRAE AAVSGDPEAP GSNTAPEPKE EGSAHLAVPG VYFTCPLTGA ILRKDQRDAR IREAILMHFS TDPVAASIMK IHTFNKDRDR VKLGVDTIAK YLDNIHLHPE EEKYRKIKVQ NKVFQERIHG LEGTHEFFEA IGFQKVLLPI PDQEGPEEFY VLSEALAQP QSLEWHKEQL LSAEPVRATL ARQRRVFRPS TLASQFDLPA DFFNLTAEEI KREQRLRSEA VERLSVLRTK AMREREEQRE MRKYTYTLRL VRLPDGCLLQ GTFYARERVA ALYGFVREAL QNDWLPEFELL ASGGQKLSL ED ENLAFNECGL VPSALLTFSL DAAVLEDIRA AGTQPDSIL KPELLSAIEK L
Source	Mammalian cell
Target Names	UBXN6
Protein Names	Recommended name: UBX domain-containing protein 6 Alternative name(s): UBX domain-containing protein 1
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.