



Recombinant Human Coiled-coil domain-containing protein 96 (CCDC96)

Product Code	CSB-MP646658HU
Abbreviation	CCDC96
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2M329
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDVSSEHTKD PGGEGGDGES LAARPSKIK SSGPPTSPEP GELESEPEEE EEEQAASQGG TAADEQAEAP KGLTAAEAAG EEGPGEPGRP AEPQPEPEEP AEVGAEPPAQ PEPGAGPEEL EAEAGAEEL QAAEGKEVRF QASLPLTRID EEEEAAPEA ETERVEGEEE DKEETQRDGA ESKERDGEGR PAKSQEEGKR LYGRDEFEDL EWSEEVQKLQ EQQLRSDLLD QYRSLVERN RSQRYNLYLQ HKIFEALRRK KGLEAAEVAD RGAEAEAPEK EQAYLRHLGM LEELKKQQAD DLQWYHQELG QLKRQCQEK TRVEKEWRRF QALKKQVVMQ AMGSCRMRRG RQAALREVEQ IQALEDKKEK EMSAVRLNI QLKQSLVHFE TRMRTQEDLT QGLLLIDFEQ LKIENQTFNE KIEERNEELL KLRSKVTNSV QVITHVKEKL HFMDMENACK KTQLAEIEAQ AALGRDILTK TKQAREGLRT DNIRLNQKCG LLGKDSLRLD LEEKVDKTEL LHRRLSLKR HHASLTLSR GVRQKIREAK AFLPS
Source	Mammalian cell
Target Names	CCDC96
Protein Names	Recommended name: Coiled-coil domain-containing protein 96
Expression Region	1-555
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.