



Recombinant Methanosphaera stadtmanae Peptide chain release factor subunit 1 (prf1)

| | |
|--------------------------|--|
| Product Code | CSB-YP646720MAAR |
| Abbreviation | prf1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q2NEL3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Methanosphaera stadtmanae (strain ATCC 43021 / DSM 3091 / JCM 11832 / MCB-3) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSDVSSKEIY EVKKTLEKELE DKKGRGTELV SVYIPPEKQI SDVAKQMRDE LGQSANIKSK QTRKNVQSAI EVIIQRLKLF PKPPEKGLVM FVGMIPKGGP GTEKMETYVF QPPEAVQTYT YHCDSQFFVE PLKQIIEYKE VYGVVVLDK ESTIATLRGK RIDIIKHLTS GVPGKHKAGG QSQRFRFDRVI ELAAHEFLKR IGRHVDEAFL PLKDELKGV L IGGPGHTKND FVDGEYIHYE IHDKIINIVD TSYTGDFGIR EVIDESADTL DEMDIMQEKK FMRKFLTGLI SESGLSTYGE KEVRQNLQMG AVETLLISEN LKSKRQTYTC PACNTVDVIT TRQHQPPEK RCPKCNEVMK ITKTQETAEE LIELAEVKT HVEVISIETE EGTQLDKAFG GIAGILRYKV K |
| Source | Yeast |
| Target Names | prf1 |
| Protein Names | Recommended name: Peptide chain release factor subunit 1 Alternative name(s): Translation termination factor aRF1 |
| Expression Region | 1-411 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.