



Recombinant Callithrix jacchus Protein Wnt-2 (WNT2)

Product Code	CSB-EP646917CYL-B
Abbreviation	WNT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QLG1
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	≥85% (SDS-PAGE)
Sequence	SWWYM GATGGSSRVM CDNVPGLVSS QRQLCHRHPD VMRAIGLGVT EWTSECQYQF RQHRWNCNTL DRDHSLFGRV LLRSSRESAF VYAISSAGVV FAVTRACSQG EVKSCSCDPK KMGSGKDSKG VFDWGGCSDN IDYGIKFARA FVDAKERK GK DARALMNLHN NRAGRKSVKR FLKQECKCHG VSGSCSLRTC WLAMADFRKT GDYLWRKYNG AIQVVMNQDG TGFTVANERF KKPTKNDLVY FENSPDYCIR DRETGSLGTA GRVCNLTSRG MDSCEVMCCG RGYDTSHVTR MIKCGCKFWH CCAVRCQDCL EALDVHTCKA PKNADWKSPT
Source	E.coli
Target Names	WNT2
Protein Names	Recommended name: Protein Wnt-2
Expression Region	26-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.