



Recombinant Human Inhibitory synaptic factor 1 (INSYN1)

Product Code	CSB-YP647179HU
Abbreviation	C15orf59
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2T9L4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNIRGAPDLG QPSDDPSSGG ERERIRQRMK MVIGQLEGIL RELKEVAKEL REVVSQIDKL TSDDFELEP DDWTTATVSS TSSSDKAGMG GPFDLGHLDLDF MTADILSDSW EFCSFLDVST PSDSVDGPES TRPGAGPDYR LMNGGTPIPN GPRVETPDSS SEEAFGAGPT VKSQLPQRTP GTRERVRFSK KVLVHALCCD DEEGDGEQEV EEEEEVGLPPE PAHTEAHAGP HKPSPAPYKS RRSPLTSRHS GSTLAPEQTR RVTRNSSTQT VSDKSTQTVL PYTATRQKAR GKN
Source	Yeast
Target Names	C15orf59
Protein Names	Recommended name: UPF0583 protein C15orf59
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.