



# Recombinant Bovine Electron transfer flavoprotein subunit beta (ETFB)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP647206BO  |
| <b>Abbreviation</b>      | ETFB  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q2TBV3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Bos taurus (Bovine)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | AELRALVAV KRVIDFAVKI RVKPKDTGVV TDGVKHSMNP FCEIAVEEAV<br>RLKEKLVKE IIAVSCGPAQ CQETIRTALA MGADRGHVE VPAAEANHLG<br>PLQVARVLAK LAEKEKVDLV LLGKQAIDDD CNQTGQMTAG FLDWPQGTF<br>SQVTLEGDKI KVEREIDGGL ETLRLKLPV VTADLRLNEP RYATLPNIMK<br>AKKKKIEVIK AGDLGVDLTS KLSVISVEDP PQRTAGVKVE TTEDLVAKLK<br>EIGRI  |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | ETFB  |
| <b>Protein Names</b>     | Recommended name: Electron transfer flavoprotein subunit beta Short name= Beta-ETF  |
| <b>Expression Region</b> | 2-255   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This gene encodes electron-transfer-flavoprotein, beta polypeptide, which shuttles electrons between primary flavoprotein dehydrogenases involved in mitochondrial fatty acid and amino acid catabolism and the membrane-bound electron transfer flavoprotein ubiquinone oxidoreductase. The gene deficiencies have been implicated in type II glutaricaciduria. Alternatively spliced transcript variants have been found for this gene. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                               |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.