



Recombinant Bovine Coiled-coil domain-containing protein 86 (CCDC86)

Product Code	CSB-EP647208BO
Abbreviation	CCDC86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TBX7
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MDTPLRRSRR LEGLKPESPE NPTSVLRVRR VLVEFESNPK ETGEPSPPPG LGSPSRQPET SPGSPSLPNG PALGSPRKQP ELDSGSPEGH RDPGLNFPQN QPESSESHL LQPKPSEESP KFSQNGQGEAD SELPKSKEEP TPGCPRHQLQ QDSGSLEFPF GQKAPGPEPS KPLQELTPRS PGSPRDQHEP SKPPAAGEPA REGPAPKKRE GSSAQAPASK KPKEEIPVIP KGKPKSGRVW KDRSKKRFSQ MVQDKPLRTS WQRKMKDRQE RKLAKDFARH LEEEKERRRQ EKKKRRAENL RRLENERKA EIVQVIRNPA KLKRAKKKQL RSIEKRDTLA QLQKQPPQRP ATKV
Source	E.coli
Target Names	CCDC86
Protein Names	Recommended name: Coiled-coil domain-containing protein 86
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.