



# Recombinant Solanum tuberosum NAD (P)H-quinone oxidoreductase subunit I, chloroplastic (ndhl)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP647288FIG   |
| <b>Abbreviation</b>      | ndhl  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q2VEC7  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Solanum tuberosum (Potato)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLPMITEFIN YGQQTIRAAR YIGQGFMITL SHANRLPVTI QYPYEKLITS<br>ERFRGRIHFE FDKCIACEVC VRVCPIDL PV VDWKLETDIR KKRLLNYSID<br>FGICIFCGNC VEYCPTNCLS MTEEYELSTY DRHELNYNQI ALGRLPMSVI<br>DDYTIRTISN LPQINNE   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | ndhl  |
| <b>Protein Names</b>     | Recommended name: NAD(P)H-quinone oxidoreductase subunit I, chloroplastic<br>EC= 1.6.5.- Alternative name(s): NAD(P)H dehydrogenase subunit I Short<br>name= NDH subunit I NADH-plastoquinone oxidoreductase subunit I  |
| <b>Expression Region</b> | 1-167   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |