



Recombinant Human Poly [ADP-ribose] polymerase 6 (PARP6)

Product Code	CSB-YP647549HU
Abbreviation	PARP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2NL67
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDIKGQFWND DDSEGDNESE EFLYGVQGSC AADLYRHPQL DADIEAVKEI YSENSVSIRE YGTIDDVDID LHINISFLDE EVSTAWKVLRLR TEPIVLRRLRF SLSQYLDGPE PSIEVFQPSN KEGFGLGLQL KKILGMFTSQ QWKHLSNDFL KTQQEKRHSW FKASGTIKKF RAGLSIFSPI PKSPSFPIIQ DSMLKGKLG PELRVGRMLN RSISCTMKNP KVEVFGYPPS PQAGLLCPQH VGLPPPARTS PLVSGHCKNI PTLEYGFLVQ IMKYAEQRIP TLNEYCVVCD EQHVFQNGSM LKPAVCTREL CVFSFYTLGV MSGAAEEVAT GAEVVDLLVA MCRAALES PR KSIIFEPYPS VVDPTDPKTL AFNPKKKNYE RLQKALDSVM SIREMTQGSY LEIKKQMDKL DPLAHPQLQW IISSNRSHIV KLPLSRLKFM HTSHQFLLLS SPPAKEARFR TAKKLYGSTF AFHGSHIENW HSILRNLGVN ASYTKLQLHG AAYGKGIYLS PISSISFGYS GMGKGQHRMP SKDELVQRYN RMNTIPQTRS IQSRFLQSRN LNCIALCEVI TSKDLQKHGN IWVCPVSDHV CTRFFFVYED GQVGDANINT QDPKIQKEIM RVIGTVQVYTN
Source	Yeast
Target Names	PARP6
Protein Names	Recommended name: Poly [ADP-ribose] polymerase 6 Short name= PARP-6 EC= 2.4.2.30
Expression Region	1-630
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.