



Recombinant Microcebus murinus Hepatocyte growth factor receptor (MET), partial

Product Code	CSB-MP647700MFU
Abbreviation	MET
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QL89
Product Type	Recombinant Protein
Immunogen Species	Microcebus murinus (Lesser mouse lemur) (Gray Mouse Lemur)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	MET
Protein Names	Recommended name: Hepatocyte growth factor receptor Short name= HGF receptor EC= 2.7.10.1 Alternative name(s): HGF/SF receptor Proto-oncogene c-Met Scatter factor receptor Short name= SF receptor Tyrosine-protein k
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.