



Recombinant *Oryza sativa* subsp. *japonica* B3 domain-containing protein Os11g0197600 (Os11g0197600, LOC_Os11g09160)

Product Code	CSB-EP647723OFG-B
Abbreviation	Os11g0197600, LOC_Os11g09160
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2R9D2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVVREKQGGR MGKKGKKGKE KECDINCFGR SFFRVLLTLQ SLERMKIPSS FNQCLQNQPT GMVSLVDRSG NKWSAELTSD SEGFFVHGW KEFVRDNSIQ CGQFLVFTYD KRSQFSVTVF EPSGIDKIST FSAHPSKNVI IKTESDEGGM VTAAITTEKM APALKENNGI TGKRTRDVDY LMEDRVVVFV KSSEANVCES SRRKRAGASA GKSKVTSTSH NSTRGSSCSS DEDNSSSKSP NPPFLMRFLS GEVSRRGRCV SKGQRQLTVI SQRRPVTEAE KDHALQRARE FKSKNPFVAVQ IMMESYVYVG FFMNIPCEV RECLPHTNKR ITLWDPQGKA WEVNYVYYSR RSVGFSFGGW GKFAVGNLLE KFDVCFVFLV QKDNKIVHIY RVVPEITPHK LRSDPK
Source	<i>E.coli</i>
Target Names	Os11g0197600
Protein Names	Recommended name: B3 domain-containing protein Os11g0197600
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.