



# Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os11g0197600 (Os11g0197600, LOC\_Os11g09160)

<b>Product Code</b>	CSB-EP647723OFG
<b>Abbreviation</b>	Os11g0197600, LOC_Os11g09160
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2R9D2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVVREKQGGR MGKKGKKGKE KECDINCFGR SFFRVLLTLQ SLERMKIPSS FNQCLQNQPT GMVSLVDRSG NKWSAELTSD SEGFFVHGW KEFVRDNSIQ CGQFLVFTYD KRSQFSVTVF EPSGIDKIST FSAHPSKNVI IKTESDEGGM VTAAITTEKM APALKENNGI TGKRTRDVDY LMEDRVVVFV KSSEANVCES SRRKRAGASA GKSKVTSTSH NSTRGSSCSS DEDNSSSKSP NPPFLMRFLS GEVSRRGRCV SKGQRQLTVI SQRRPVTEAE KDHALQRARE FKSKNPFVAVQ IMMESYVYVG FFMNIPCEV RECLPHTNKR ITLWDPQGKA WEVNYVYYSR RSVGFSFGGW GKFAVGNLLE KFDVCFVFLV QKDNKIVHIY RVVPEITPHK LRSDPK
<b>Source</b>	E.coli
<b>Target Names</b>	Os11g0197600
<b>Protein Names</b>	Recommended name: B3 domain-containing protein Os11g0197600
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.