



Recombinant Human Peroxisomal leader peptide-processing protease (TYSND1)

Product Code	CSB-MP647967HU
Abbreviation	TYSND1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2T9J0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRRQWGSAMR AAEQAGCMVS ASRAGQPEAG PWSCSGVILS RSPGLVLCHG GIFVPFLRAG SEVLTAAGAV FLPGDSCRDD LRLHVQWAPT AAGPGGGAER GRPGLCTPQC ASLEPGPPAP SRGRPLQPRL PAELLLLSC PAFWAHFARL FGDEAAEQWR FSSAARDDEV SEDEEADQLR ALGWFallGV RLGQEEVEEE RGPAMAVSPL GAVPKGAPLL VCGSPFGAFC PDIFLNTLSC GVLSNVAGPL LLTDARCLPG TEGGGVFTAR PAGALVALVV APLCWKAGEW VGFTLLCAAA PLFRAARDAL HRLPHSTAAL AALLPPEVGV PWGLPLRDSG PLWAAAVALV ECGTVWGSV AVAPRLVVC RHVSPREAR VLVRSTTPKS VAIWGRVVFA TQETCPYDIA VVSLEEDLDD VPIPVPAEHF HEGEAVSVVG FGVFGQSCGP SVTSGILSAV VQVNGTPVML QTTCVHSGS SGGPLFSNHS GNLLGIITSN TRDNNTGATY PHLNFSIPIT VLQPALQQYS QTQDLGGLRE LDRAAEPVRV VWRLQRPLAE APRSKL
Source	Mammalian cell
Target Names	TYSND1
Protein Names	Recommended name: Peroxisomal leader peptide-processing protease EC=3.4.21.- Alternative name(s): Trypsin domain-containing protein 1 Cleaved into the following 2 chains: 1. Peroxisomal leader peptide-processing protease, 15 kDa form
Expression Region	1-566
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.