



Recombinant Arabidopsis thaliana F-box protein At3g59000 (At3g59000)

Product Code	CSB-BP648067DOA
Abbreviation	At3g59000
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2V3N5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDRVGSLPDE LLSHILSFLT TKEAALTSLL SKRWRYLIAF VPNLAFDDIV FLHPEEGKPE RDEIRQSFMD FVDRVLALQA ESPIKKFSLK CRIGVDSDRV DGWISNVLNR GVSELDLLII LGMTMEDSYR LSPKGFASKT LVKLEIGCGI DISWVAGSIF LPMLKTLVLD SVWFYVDKFE TLLLALPALE ELVLVDVNW DSDVTISNAS LKTLTIDSDG HLGTFSDTP SLVYFCYSYD AAEDYPVVKM ENLREARIFL LVTDNEIERV RLPINDGLED AMDNVVLRFG HVGKLMNGIR NVEYLELSAD TLEVLSLCCE SMPVFKNLKS LAIKSTEGRG WQAMPVLLRN CPHLEFLIE GLLHHVTDKC GDACDCVPRE DKGRSLTSCP VNMLEIHGFR GTTKEMQMIK HFLDYFPSLK EMDIYADEDG PTNLEAPGMF EQITQLFTLY DEFYTCDVQF MVRGSLYKKL TAQCCVFKEE NLPHRLVRTL G
Source	Baculovirus
Target Names	At3g59000
Protein Names	Recommended name: F-box protein At3g59000
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.