



Recombinant *Oryza sativa* subsp. japonica Putative germin-like protein 12-3 (Os12g0154900, LOC_Os12g05870)

| | |
|--------------------------|---|
| Product Code | CSB-EP648512OFG |
| Abbreviation | Os12g0154900, LOC_Os12g05870 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q2QXJ1 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SDPSPLQD LCVADKNSPV RVNGFPCKDA KDVSVDFFL AANLDKPMDI TKSKAGSNVT LINVMKLAGL NTLGISMARI DYAPKQNP HTHPRATEIL SVIEGSLYVG FVTSNQANGE NKLFTKTLNK GDVFVPEGL IHFQFNPSYD KPAAAIVALS SQNPGAITIA NAVFGSNPPI SDDVLAKAFQ VDKKAVDWLQ AQFWENNHN |
| Source | E.coli |
| Target Names | Os12g0154900 |
| Protein Names | Recommended name: Putative germin-like protein 12-3 |
| Expression Region | 23-229 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |