



Recombinant *Oryza sativa* subsp. *japonica* Mannan endo-1,4-beta-mannosidase 7 (MAN7)

Product Code	CSB-BP648521OFG
Abbreviation	MAN7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2RBB1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEDAHHHHWT MVERRGSQLW ASGRPFIHG FNTYWLMSFA ADQATRLRVT AAIAEAGLNV CCTWAFSDGG YRALQTAPFH YDEDVFRALD FVNSEARRHN MRLILSLCNN WEDYGGKAQY VRWGKEAGLD LTSEDDFFSD PTIKSYYKAF VEAVVTRINT VTNETYKDDP TILAWELINE PRCPSPSGD TLQAWMEEMA SYVKSIDPVH LLEIGIEGFY GPSIPPELLPV NPDEYSGHAG IDFIRNHQAP GIDLASIHVY SDIWL PQSIK ENHLQFVDKW MQQHIDDAAN LLGMPIVVGE FGVSVKDGKF GNEFREDFMK TIYRIFLSSW KEGVIGGGCL LWQLFPEGAE HMDDGYAVIF AKSPSTLSLL ANHLRCLEC
Source	Baculovirus
Target Names	MAN7
Protein Names	Recommended name: Mannan endo-1,4-beta-mannosidase 7 EC= 3.2.1.78 Alternative name(s): Beta-mannanase 7 Endo-beta-1,4-mannanase 7 OsMAN7
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.