



Recombinant Mouse Protein C-ets-2 (Ets2)

Product Code	CSB-MP007853MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15037
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNDFGIKNMD QVAPVANSFR GTLKRQPAFD TFDGSLFAVL PSLSEDQTLQ EVPTGLDSVS HDSASCELPL LTPCSKAVMS QALKATFSGF QKEQRRLGIP KNPWLWSEQQ VCQWLLWATN EFSLVNVNLH QFGMNGQMLC NLGKERFLEL APDFVGDILW EHLEQMIKEN QEKTEDQYEE NSHLNAVPHW INSNTLGFMS EQAPYGMQAP NYPKDNLLDS MCPPTSATPAA LGSELQMLPK SRLNTVNVNY CSISQDFPSS NVNLLNNSG KPKDHDSPEN GGDSFESSDS LLRSWNSQSS LLDVQRVPSF ESFEEDCSQS LCLSKLTMSF KDYIQERSDP VEQGKPVIPA AVLAGFTGSG PIQLWQFLE LLSDKSCQSF ISWTGDGWEF KLADPDEVAR RWGKRKNKPK MNYEKLRSGL RYYYDKNIIH KTSKGKRYVYR FVCDLQNLLG FTPEELHAIL GVQPDTED
Source	Mammalian cell
Target Names	Ets2
Protein Names	Recommended name: Protein C-ets-2
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.