



Recombinant Human PWWP domain-containing protein MUM1 (MUM1)

Product Code	CSB-MP648777HU
Abbreviation	MUM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TAK8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MADAKYVLCR WEKRLWPAKV LARTATSTKN KRRKEYFLAV QILSLEEKIK VKSTEVEILE KSQIEAIASS LASQNEVPAA PLEELAYRRS LRVALDVLSE GSIWSQESSA GTGRADRSR GKPMEHVSSP CDSNSSSLPR GDVLGSSSRPH RRRPCVQQSL SSSFTCEKDP ECKVDHKKGL RKSENPRGPL VLPAGGGAQD ESGSRIHHKN WTLASKRGGN SAQKASLCLN GSSLSEDDTE RDMGSKGGSW AAPSLPSGVR EDDPCANAEG HDPGLPLGSL TAPPAPEPSA CSEPGCEPAK KRPRLDGSQR PPAVQLEPMA AGAAPSPGPG PGPRESVTPR STARLGPPPS HASADATRCL PCPDSQKLEK ECQSSEESMG SNSMRSILEE DEEDEEPPRV LLYHEPRSFE VGMLVWHKHK KYPFWPAVVK SVRQRDKKAS VLYIEGHMNP KMKGFTVSLK SLKHFDCKEK QTLLNQARED FNQDIGWCVS LITDYRVRLG CGSFAGSFLE YYAADISYPV RKSIIQQDVLG TKLPQLSKGS PEEPVVGCP L GQRQPCRKML PDRSRAARDR ANQKLVEYIV KAKGAESHLR AILKSRKPSR WLQTFLSSSQ YVTCVETYLE DEGQLDLVVK YLQGVYQEVG AKVLQRTNGD RIRFILDVLL PEAIICAISA VDEVYKTAE EKYIKGPSLS YREKEIFDNQ LLEERNRRRR
Source	Mammalian cell
Target Names	MUM1
Protein Names	Recommended name: PWWP domain-containing protein MUM1 Alternative name(s): Mutated melanoma-associated antigen 1 Short name= MUM-1 Protein expandere
Expression Region	1-710
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.