



# Recombinant *Aspergillus oryzae* Autophagy protein 5 (atg5)

<b>Product Code</b>	CSB-EP648816DPA-B
<b>Abbreviation</b>	atg5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2UBM1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Aspergillus oryzae</i> (strain ATCC 42149 / RIB 40) (Yellow koji mold)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METQATLNSI QKAVWDGRLP LQIRLAPSES RIYDQTDPYL ISYPRISYLP SLLPRLRAFF ASSLIDPSSN AHDGWFSFEG VPLKWHLPVIG LLYDLYAGAD PASKGTAESE DAGWDIDDQD NPLPWRLVVH FSDWPDEELV RLDAEGMVMN DAFINSVKEA DFLRNGTAKG IMSLSKEDSS GLWKSVQNVE LSSFQRISNI LLPPLNQPF R NIPRIFLPL PPDSGSPSLK IVQSPVPLI PPSSVAASQL ALSRSSITPQ TQTIGSALHS LLPNLFPSRR TPVLAKPVLH GAAVPMSAPV EELVRSSAYG DGWLYVVIRM MG
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	atg5
<b>Protein Names</b>	Recommended name: Autophagy protein 5
<b>Expression Region</b>	1-322
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.