



Recombinant Human Serine/threonine-protein kinase B-raf (BRAF)

Product Code	CSB-EP002791HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15056
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>AALSGGGGG GAEPGQALFN GDMEPEAGAG AGAAASSAAD PAIPEEVWNI KQMIKLTQEH IEALLDKFGG EHNPPSIYLE AYEEYTSKLD ALQQREQQLL ESLGNGTDFS VSSASMDTV TSSSSSSLSV LPSSLSVFQN PTDVARSNPK SPQKPIVRVF LPNKQRTVVP ARCGVTVRDS LKKALMMRGL IPECCAVYRI QDGEKKPIGW DTDISWLTGE ELHVEVLENV PLTTHNFVRK TFFTlafCDF CRKLLFQGFR CQTCGYKFHQ RCSTEVPLMC VNYDQLDLLF VSKFFEHHPI PQEEASLAET ALTSGSSPSA PASDSIGPQI LTSPSPSKSI PIPQFRRPAD EDHRNQFGQR DRSSAPNVH INTIEPVNID DLIRDQGFGRG DGGSTTGLSA TPPASLPGSL TNVKALQKSP GPQRERKSSS SSEDNRNRMKT LGRRDSSDDW EIPDGQITVG QRIGSGSFGT VYKGKWHGDV AVKMLNVTAP TPQQLQAFKN EVGVLKTRH VNILLFMGYS TKPQLAIVTQ WCEGSSLYHH LHIIETKFEM IKLIDIARQT AQGMDYLHAK SIIHRDLKSN NIFLHEDLTV KIGDFGLATV KSRWSGSHQF EQLSGSILWM APEVIRMQDK NPYSFQSDVY AFGIVLYELM TGQLPYSNIN NRDQIIFMVG RGYLSPDLK VRSNCPKAMK RLMAECLKKK RDERPLFPQI LASIELLARS LPKIHRASE PSLNRAGFQT EDFSLYACAS PKTPIQAGGY GAFPVH</p>
Source	E.coli
Target Names	BRAF
Protein Names	Recommended name: Serine/threonine-protein kinase B-raf EC= 2.7.11.1 Alternative name(s): Proto-oncogene B-Raf p94 v-Raf murine sarcoma viral oncogene homolog B1
Expression Region	2-766
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein belonging to the raf/mil family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene are associated with cardiofaciocutaneous syndrome, a disease characterized by heart defects, mental retardation and a distinctive facial appearance. Mutations in this gene have also been associated with



various cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, thyroid carcinoma, non-small cell lung carcinoma, and adenocarcinoma of lung. A pseudogene, which is located on chromosome X, has been identified for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.