



Recombinant Horse Protein Wnt-2 (WNT2)

Product Code	CSB-YP649285HO
Abbreviation	WNT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QLA5
Product Type	Recombinant Protein
Immunogen Species	Equus caballus (Horse)
Purity	>85% (SDS-PAGE)
Sequence	SWWYM RATGGASRVM CDNVPGLVSR QRQLCHRHPD VMRAIGLGVA EWTAECQHQQF RQHRWNCNTL DRDHSLFGRV LLRSSRESAF VYAISSAGVV FAITRACSQG ELRSCSCDPK KKGTAKDSKG TFDWGGCSDN IDYGIKFARA FVDAKERKGG DARALMNLHN NRAGRKAVKR FLKQECKCHG VSGSCTLRTC WLAMADFRKT GDYLWRKYNG AIQVVMNQDG TGFTVANKKF KKPTKNDLVY FENSPDYCIR DRDAGSLGTA GRVCNLTSRG MDSCEVMCCG RGYDTSRVTR MTKCECKFW CCAVRCQDCL EALDVHTCKA PKSADWASPT
Source	Yeast
Target Names	WNT2
Protein Names	Recommended name: Protein Wnt-2
Expression Region	26-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.