



Recombinant *Oryza sativa* subsp. japonica Probable calcium-binding protein CML25/26 (CML25)

Product Code	CSB-BP649301OFG
Abbreviation	CML25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QYW1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASSASSVFA AFDKDGDKV SASELRGCM AALGEEVSEE EAAAILATAD TDGDGLLDHH EFMRLSAAHQ LQEPAEESLR CLREAFDMYA EEEETAVITP ASLRRMLRRL GSEHQRLME ECRAMICRFD LNGDGVLSFD EFRVMMLMA
Source	Baculovirus
Target Names	CML25
Protein Names	Recommended name: Probable calcium-binding protein CML25/26 Alternative name(s): Calmodulin-like protein 25/26
Expression Region	1-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.