



Recombinant *Oryza sativa* subsp. japonica Germin-like protein 11-1 (Os11g0537300, LOC_Os11g33110)

Product Code	CSB-YP649305OFG
Abbreviation	Os11g0537300, LOC_Os11g33110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2R352
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	DSPPLQD VCPMAPQGER KLFMNGFFCK SPSTIMASDF KTLNHNHAGD LDNMVRRSSAN IITATEFPGL NTLGISMART DIAVSGAVLP HSHPRASEMM FVHSGSVVAG FFDTKGKLFQ KTLAEGDVFI FPRGLVHFIM NYGFGLATTF SVLNSQNPV VGITHAMFAP DSEVAEGLMA RMLSFDRDGM DDSSVDSPW FY
Source	Yeast
Target Names	Os11g0537350
Protein Names	Recommended name: Germin-like protein 11-1
Expression Region	24-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.