



Recombinant Human Ezrin (EZR)

| | |
|--------------------------|---|
| Product Code | CSB-EP007914HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P15311 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | PKPINVRVT TMDAELEFAI QPNTTGKQLF DQVVKTIGLR EVWYFGLHYV DNKGFPTWLK LDKKVSAQEV RKENPLQFKF RAKFY PEDVA EELIQDITQK LFFLQVKEGI LSDEIYCPPE TAVLLGSYAV QAKFGDYNKE VHKSGYLSSE RLIPQRVMDQ HKLTRDQWED RIQVWHA EHR GMLKDNAMLE YLKIAQDLEM YGINYFEIKN KKGTDLWLGV DALGLNIYEK DDKLTPKIGF PWSEIRNISF NDKKFVIKPI DKKAPDFVFY APRLRINKRI LQLCMGNHEL YMRRRKPDTI EVQQMKAQAR EEKHQKQLER QQLETEKKRR ETVEREKEQM MREKEELMLR LQDYEEKTKK AERELSEIQ RALQLEEERK RAQEEAERLE ADRMAALRAK EELERQAVDQ IKSQEQLAAE LAEY TAKIAL LEEARRRKED EVEEWQHRAK EAQDDL VTK EELHLVMTAP PPPPPVYEP VSYHVQESLQ DEGA EPTGYS AELSSEGIRD DRNEEKRITE AEKNERVQRQ LLTLSS ELSQ ARDENKRTHN DIIHNENMRQ GRDKYKTLRQ IRQGNTKQRI DEFEAL |
| Source | E.coli |
| Target Names | EZR |
| Protein Names | Recommended name: Ezrin Alternative name(s): Cytovillin Villin-2 p81 |
| Expression Region | 2-586 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.